

Portable Colorimetric Analyser

PRODUCT DATASHEET

APPLICATIONS

Wastewater Process Water Drinking Water Surface Water Seawater

MEASUREMENTS

Aluminium Ammonia Iror Manganese Nicke Nitrate Nitrite Orthophosphate TON

FEATURES

Flexible Loop Flow Analysis (LFA)*

Multi-parameter Options

*Patented by Systea Srl, Italy

INSTALLATION OPTIONS

Work Top
Field Deployable Weatherproof
Package
Power options – mains or battery with
renewable energy options








The MicroMac 1000 is designed to operate as a free standing colorimetric analyser which gives the user the ability to get "live" data from a site for a short period of time without having to install complex on-line monitoring packages.

The portability of the analyser gives the user a portable "laboratory" as it uses the same standard methods methods of analysis as most laboratories. It can be used as a tool for site surveys, investigations or when set up as a package, longer term studies or profiling of a site. With the advantage of "live" results the system can be used in place of samplers where the sample is collected and sent away for laboratory analysis thus enabling changes to be made much quicker.

The chemistry module employs the patented Loop Flow Analysis System (LFA –Systea Srl) resulting in a very flexible chemistry base which can incorporate high temperature digestion or heating, UV digestion, acid digestion, temperature controlled end point reactions all utilising either a visible light colorimeter or a fluorimeter.

In the same way as the wall mounted monitor, the MicroMac C, the MicroMac 1000 can be configured as a single or multiple chemistry system up to a maximum of 4





Call us on 01726 879800 www.partech.co.uk






Loop Flow Analysis

PRODUCT DATASHEET

The analytical methods used in the MicroMac 1000 are based on the Standard Methods that are used in Laboratories throughout the world. Many of these methods have their basis in the "Methods for the Examination of Waters and Associated Materials" as published by HMSO – Standing Committee of Analysts.

The Loop Flow Analysis system can be configured as a single or multiple chemistry module: The multiple chemistry module can have up to 4 chemistries or 7 reagent additions. The multiple chemistry module above shows the inclusion of a high temperature heating bath for hot acid digestion and has high pressure valves and teflon tubing in all areas in contact with the acid or high temperatures.





STANDARD METHODS (others available on request)

Aluminium 0-500 μg/l to 0-5 mg/l Orthophosphate 0-500 μg/l as P to 0-50 mg/l as P

Ammonia 0-2 to 0-50 mg/l as N Soluble Iron 0-500 μg/l to 0-10 mg/l

Manganese 0-300 μg/l to 0-5 mg/l

Nitrate (TON) 0-1 to 0-20 mg/l as N

Nitrite 0-1 to 0-20 mg/l as N

Call us on 01726 879800 www.partech.co.uk



Field Deployable Package


PRODUCT DATASHEET



This site evaluation package allows the user to quickly assess problems that a site is encountering by monitoring 24 hours per day and at weekends. Data from site evaluations has already allowed water company's to identify site process problems and to track down illegal dumping of chemicals into the water course by an unscrupulous manufacturer.


Partech are offering this package for hire in a package that includes the analyser, sample pump, delivery and setup, reagents and data reporting. This package can be tailored to suit the particular needs of a site.

Alternatively the MicroMac I 000 Portable analyser can purchased along with its sister product, the MicroMac C which is designed for permanent installation.



Result Set showing two large Phosphate peaks which were subsequently traced to a Trade Effluent discharge









Portable - Colorimetric Analyser

PRODUCT DATASHEET

Analyser

Measuring Principle

Colorimeter

Fluorimeter

Number of Parameters

Measurement Frequency

Measurement Time

Number of Sample Points

Sample Requirements

Waste

Reagent Cooler

Physical

Mounting

Protection Rating

Weight

Dimensions

Environmental Temperature

Electrical

Power Supply

Power Use

Hardware

Communication Port

Output Signals

Input Signals

Alarm Signals

Alarm Messages

Colorimetric or Fluorimetric Dual Beam, Silicon Detector

Excitation at 370 nm, emission 420-470 nm

I standard, upto 4 depending on combination

Programmable

Method Specific

I standard, upto 6 optional

10 to 30 C

Toxic and Non Toxic fed to separate drain

optional Pelter Cell

Wall Mounting normally in building or kiosk

IP55

25 kg without reagents

 $800 \times 450 \times 300 \text{ mm (hxwxd)}$

10 to 30 C

12VDC or 115/230 VAC

Typically 4W on standby, 10W during analysis

PC104 industrial standard, integrated keyboard and display

RS232, RS485, USB

4-20mA per parameter, 400 ohm maximum load

Remote analysis and calibration request

1x High Alarm, SPDT, 24VDC, 0.5A per parameter

 $I \times General Alarm, SPDT, 24VDC, 0.5A$

Ix Calibration Alarm, SPDT, 24VDC, 0.5 A per parameter

On display

Publication No: 184640DS-Iss02 The company reserves the right to alter the specification without prior notice. E&OE

